


MagicPure™ Viral DNA/RNA Kit

MagicPure™ Viral DNA/RNA Kit is designed for DNA/RNA isolation from up to 200 µl of plasma, serum, body fluid and mammalian cell culture supernatant samples. A workflow based magnetic-bead technology guarantees a high-throughput scale and speed with yield and purity. Isolated DNA/RNA is free of protein contamination, and it is suitable for a broad range of applications, including PCR, RT-PCR, qPCR and sequencing.

- No centrifugation required, suitable for both manual and automated DNA/RNA isolation.
- Fast and easy, with enhanced purity and yield.



RNA isolated from BVDV, and then detected by qRT-PCR
Green: TransGen
Red: Company TH